



Microsynth

Overnight Sequencing Service



For Barcode Economy Run and Economy Run samples.



Find the mailbox nearest to you: Login to our Homepage (-> View Collection Points).



Take advantage of our FREE collection service: we bring your samples to Microsynth at no additional cost, using environmentally friendly transport (bicycle/train).



You receive the sequencing results in the morning the following business day.

Instructions for Barcode Economy or Economy Run samples

Individual Sequences with reads up to 1000 bases

1 tube per reaction, total volume: 10 µl

- 0.8 µg Plasmid DNA or 15ng/100 bases PCR Product**
- 20 pmol specific sequencing primer (standard primers are added by us)**

Accreditation according to EN ISO 17025
Certification according to ISO 9001:2000